

EUROFIX

CRANIOPLASTIE / CRANIOPLASTY



EUROFIX CRANIO

Ciment chirurgical stérile pour la cranioplastie et reconstructions osseuses.
Surgical cement sterile for cranioplasty and bone repairs.

Spécialement conçu pour être utilisé en:

- Neurochirurgie.
- Chirurgie plastique.
- Chirurgie maxillo-faciale.
- Oto-rhino-laryngologie.

Principaux avantages par rapport aux ciments utilisés en orthopédie:

- Particules de PMMA plus petites.
- Empreintes plus faciles à réaliser.
- Temps de manipulation et de modelage améliorés.
- Fini de surface plus lisse.

Méthodes de fixation:

- Matériel de suture non résorbable.
- Fils ou câbles en titane.
- Systèmes de plaque / vis en titane.

Présentation:

EUROFIX Cranio est conditionné dans un emballage contenant deux doses de ciment chirurgical en blisters indépendants; chacune de ces doses comportant: un sachet avec 30 g de poudre stérilisée à l'oxyde d'éthylène et une ampoule avec 17 ml de liquide stérilisé par ultrafiltration.

LE CIMENT EUROFIX CRANIO EST UTILISÉ DANS LA RÉPARATION DES DÉFAUTS OSSEUX D'ORIGINE TRAUMATIQUE, TUMORALE OU CHIRURGICALE.

Specially formulated for:

- Neurosurgery.
- Plastic surgery.
- Maxillo-facial surgery.
- ENT surgery.

Main advantages in comparison to orthopedic cements:

- Smaller PMMA particles.
- Improved manipulation and modeling time.
- Better impression taking.
- Smoother finishing.

Fixation methods:

- Non reabsorbable surgical sutures.
- Titanium wires or cables.
- Titanium plates / screws.

Presentation:

EUROFIX Cranio is supplied in a container with two doses of surgical cement in separate blisters, made up of: a pouch with 30 g of powder sterilized by ethylene oxide and a blister with 17 ml of liquid sterilized by ultrafiltration.

EUROFIX CRANIO CEMENT IS INTENDED FOR THE REPAIR OF TRAUMATIC, TUMORAL OR SURGICAL BONE DEFECTS.



REF 880135



MEDICAL DEVICES DESIGNER AND
MANUFACTURING SPECIALIST.
www.synimed.com

DISTRIBUTOR: _____



Synergie Ingénierie Médicale S.A.R.L.

Zone Artisanale de l'Angle - 19370 Chamberet - France - Tél.: (+) 33 5 55 98 31 38 / Fax: (+) 33 5 55 98 35 76

970664/F5/01N